

All data taken at Pacific Northwest National Laboratory (PNNL)

Operators: Steven W. Sharpe, Timothy J. Johnson and Robert L. Sams : [sw.sharpe@pnl.gov](mailto:sw.sharpe@pnl.gov)

Version 1.0, October, 01

Composite spectrum for EDB\_25T

Effective burden of composite spectrum: 1 part-per-million-meter (ppm-meter) at 296 K

Equivalent concentration x path-length of composite spectrum:  $7.734 \times 10^{-6}$  grams/liter-meter

### Sample Conditions-

- Chemical name and CAS number: 1,2-Dibromoethane, EDB, dibromoethane, Ethylene Bromide, Glycol Dibromide, sym-dibromoethane, Dowfume W 85, Bromofume, Celmid, E-D-Bee, Kopfume, Nephis, alpha,beta-dibromoethane, glycol bromide, DBE, dowfume 40, dowfume edb, dowfume w-8, dowfume w-90, dowfume w-100, EDB-85, fumo-gas, iscobrome d, pestmaster, pestmaster edb-85, soilbrom-40, soilbrom-85, soilbrom-90, soilbrom-100, soilfume, unifume, Aadibroom, Dibromoethylene, Garden dowfume, Soilbrom 90EC, Dibromoethane, 1,2- , 1,2-Dibromomethane,  $\text{CH}_2\text{BrCH}_2\text{Br}$ : [106-93-4]
- Physical properties: fw=187.862 g/mole, fp=9.97° C, bp=131.7° C
- Supplier and stated purity: Aldrich, 99+%
- Sample class: I (PNNL scale).
- Temperature of sample:  $25.06 \pm 0.02$  C
- Diluent: Sample back filled with ultra high purity nitrogen to  $760 \pm 5$  Torr
- Individual samples at 1.2090, 5.4134, 0.58360, 2.0487, 3.0645, 4.3138, 10.1400, 7.1610, 8.6501 and 6.4920 Torr. Path length = 19.96 cm. Final data is a composite spectrum.
- Preparation: Multiple freeze-thaw cycles at 0 C to remove air.

### Instrument Parameters-

- Bruker-66V FTIR, temperature controlled environment, evacuated optics bench
- Modified to include second aperture, between interferometer output and sample cell. This substantially reduces both “ghosting” and warm aperture effects.
- Spectral range: 6,500 to 570  $\text{cm}^{-1}$  (1.534 to 17.544 microns)
- Instrumental resolution based on maximum interferometer displacement is 0.112  $\text{cm}^{-1}$
- Spectral interval after 2X zero-filling interferogram and FFT: 0.06  $\text{cm}^{-1}$
- Interferogram zero-fill: 2X
- Apodization: Boxcar
- Phase correction: Mertz
- Beam splitter: Potassium bromide (KBr)
- IR source: Carbide glowbar (22 V)
- Scanner velocity: 60KHz (HeNe crossing frequency)
- Number of interferograms averaged per single channel spectra: 256
- Detector: Mid-band HgCdTe, photoconductive, 77K operation
- Folding limits: 15798 to 0  $\text{cm}^{-1}$

### Post Processing and Related Parameters-

- Non-linearity detector correction (Bruker proprietary) applied to interferogram ( $\alpha = 0.90$ ,  $\beta = 500$ )
- Composite spectrum created from 10 individual absorbance (base-10) spectra via classical least squares fit: Intercept=0, slope is fitted, individual absorbance values weighted by  $T^2$  (transmission squared), all absorbance values  $> 1.6$  are given zero weight
- Calculated and estimated errors: Type A = 1.07%, Type B = 3%

- Frequency correction (already applied):  $V(\text{corrected}) = V(\text{instrument}) * 0.999998 + 1.287 \times 10^{-4}$
- Axis units: X=wavenumbers ( $\text{cm}^{-1}$ ), Y=Absorbance (base-10)
- Trace water vapor features removed via spectral subtraction
- Baseline correction via 7<sup>th</sup> order polynomial subtraction